



**ATTENDANT
CONTROL.**

INDEPENDENCE. WITH SOME ASSISTANCE: CO-PILOT 2.

Attendant control is the extra help that is sometimes needed to enable you to participate in social activities and other aspects of an active life. The Co-pilot 2 is designed to be used often, by various people in your daily life: Family and friends, personal assistants and other care-givers, transport personnel and drivers, teachers and many others.

That's why we've made it easy to start, quick to adjust and convenient to use. The wheelchair is easy to maneuver since it adjusts to the attendants walking speed and can "feel" every change of direction through sensors in the handle. The Co-pilot 2 is used as an attendant control when needed and can be folded away neatly when not in use. Alternating

between the normal joystick and the attendant control can be done from the control panel.

It is designed for use with Corpus® seats that have R-net controls.

FEATURES



Intuitive driving.

Co-pilot 2 gives the attendant full control of the wheelchair. The handle senses your movements and adjusts the speed and direction of the wheelchair.



Safety first.

The handle functions are automatically turned off when folded. The slow speed mode makes it easier to maneuver in tight environments.



Ergonomic.

The specially shaped soft rubber handle, together with the low pressure safety grip and easy adjustment, gives an ergonomic and comfortable driving position for every user. The safety grip is designed to function with gloves in cold weather.



Simple controls.

Control panel shows a clear view of the wheelchair's status and functions: battery level, on/off, slow indoor mode, lights and alarm indicators.

BENEFITS

Easy to use

The Co-pilot 2 is intuitive to get started. Just power on, press the safety grip and slightly push forward and off you go! It requires hardly any force to control the wheelchair. The sensors in the handle operate the wheelchair motors, so you can steer the chair smoothly and predictably.

Reduce fatigue

With the smart and quick adjustments it is easy to find the best position of the handle that fits each attendant. Adjustments are important for the attendant to find the best ergonomic position for long walks, to help prevent tired arms and shoulders. The ergonomic handle with the smart safety grip also reduces the risk of fatigue.



permobil