

THE ROAD TO INDEPENDENCE.


CONTINUOUSLY STRIVING TO ENHANCE THE MOBILITY OF
WHEELCHAIR USERS, PERMOBIL PROUDLY INTRODUCES PERMOLOCK C

PERMOLOCK C



The Art of Mobility

PERMOLOCK C ALLOWS USERS TO DRIVE WHILE SITTING IN THEIR WHEELCHAIRS.



RETRACTABLE LOCKING PIN
When activated by the Control Button, the retractable locking pin drops into place to secure the wheelchair to the base. When not in use, the pin retracts for maximum ground clearance.



C400 Corpus



C500 Corpus

PERMOLOCK C
Permlock C is designed only for use with the C400 and C500 power wheelchairs.



CONTROL BUTTON
Located on the armrest, the control button is used to drop the locking pin into place, so it can properly dock to the locking plate.



LOCKING PLATE
An easy-to-install docking station that holds the chair firmly and securely in position.



VEHICLE CONTROL BUTTON
Located on the dash of the vehicle, the control button uses three lights—red, amber and green—to tell the user whether or not the chair is locked into place. To unlock the wheelchair from the locking plate, simply press and hold the control button until the locking mechanism is released.



EMERGENCY RELEASE HANDLE
Locking mechanism can be opened manually with the emergency handle on the front edge of the locking plate.



MANUAL RELEASE
Wheelchair can be released manually in the event of vehicle power failure.

ACCESSORIES



CHAIR BASE
Permits use of original car seat in Permlock C.




QUICK ANCHORING SET
Used to anchor the Permlock C in a vehicle and permits easy transfer of the lock between different positions in the vehicle.



INSTALLATION SPACER
Used for installing Permlock C on a sloping surface.

PERMOLOCK C



**PERMOLOCK C IS DESIGNED ONLY FOR USE
WITH THE C400 AND C500 POWER WHEELCHAIRS.**



The Art of Mobility

Permobil Ltd.
Unit 4, West Vale Building
Wakefield Road
Brighouse HD6 1PE
Tel 01484 722 888
Fax 01484 723013
zahida@permobiluk.co.uk

Permobil Europe BV
De Doom 22
6419 CX Heerlen
Nederland
Tel +31 45 564 54 90
Fax +31 45 564 54 91
europe@permobil.nl

www.permobil.com