

SmartDrive

Featuring
SwitchControl

permobil





SwitchControl

With the push of a button, SmartDrive can be engaged for a momentary burst of power or activated in latched mode for consistent power over extended distances. It also works as a shut off switch while the motor is running. When touched, the motor will disengage and the user stops by grabbing the wheels.

SmartDrive Specifications

Dimensions	Length: 389 mm Width: 141 mm Height: 242 mm
Product weight	5.7 kg Drive Unit
User weight	14 to 150 kg
Driving range*	Up to 19.8 km Level/No Resistance
Motor power	250 W Brushless DC
Maximum speed	6 km/h
Operating temp	-25°C to 50°C
Chair type compatibility	Rigid, Folding, Tilt-In Space, One-Arm Drive and others
Drive wheel diameter	501, 540, 559, 590 mm
Warranty	2-Year Limited Warranty

Bluetooth® LE	FCC ID IC ID	QOQBT113 5123A-BGTBLE113
---------------	-----------------	-----------------------------

Local reimbursement	Code
---------------------	------

**Performance measurements will vary based on user and chair characteristics, driving and battery conditions.*

The SmartDrive has been tested and conforms to all applicable requirements of ANSI/RESNA Standards for Wheelchairs - Volume 1 and Volume 2 and EN 12184.

NOTE: Wheelchair specifications provided by its manufacturer may be slightly affected by the addition of the SmartDrive.

Please refer to your owner's manual for more detailed operation and care instructions.