

# SmartDrive

Featuring the  
PushTracker E2



permobil





Touchscreen technology makes the SmartDrive as flexible, accessible and convenient as possible

## PushTracker E2

SmartDrive gives manual wheelchair users more control over their power assist experience than ever before. The PushTracker E2 makes using your SmartDrive even easier, with gesture-activated controls and activity tracking technology.

8x more reliable Bluetooth® connection

## Specifications

Dimensions	46.9 x 52.2 x 12.9 mm
Colors	Black
Operating system	Wear OS by Google™
Phone compatibility	Android™, iPhone™
Nominal operating voltage	3.9 V
Electric charge	415 mA-h
Wireless operating band	Bluetooth / BLE 2.4 GHz (2.402 to 2.480 GHz) WiFi 802.11 b/g/n 2.4 GHz (2.412 to 2.472 GHz) GPS / GLONASS / Beidou 1575 MHz (1.559 to 1.610 GHz)
Contains Bluetooth® LE	FCC ID 2AP42-WG12026 IC ID 24006-WG12026
Water resistance	5 ATM (Swim-ready, up to 50 m)
NFC payment	No

Band width	22 mm
Band material	Silicone (interchangeable) - Latex-Free and completely Biocompatible
Case material	Polycarbonate
Display	AMOLED (400 x 400 pixels)
App version(s)	Check Google Play™ Store

*All figures shown are for illustration purposes only and actual products may vary.  
Wear OS by Google, iPhone, Android, Bluetooth and Google Play are trademarks  
or registered trademarks of their respective owners.*