

permobil

SELECT SERIES
CUSHIONS

ROHO[®]



**EXTRAORDINARY.
JUST LIKE YOU.**

SUPERIOR SKIN & SOFT TISSUE PROTECTION HAS NEVER FELT BETTER

The SELECT Series from ROHO®

The SELECT Series is one of ROHO's many popular product lines, with the revolutionary ISOFLO® Memory Control which offers advanced skin and soft tissue protection as well as positioning capabilities. The ISOFLO provides each individual an adjustable, quick and simple stability and positioning solution while seated.

BENEFITS

Skin & Soft Tissue Protection

Based on over 45 years of research and individual experience, ROHO® DRY FLOATATION® Technology provides unparalleled performance in facilitating blood flow by fitting and matching the shape of an individual.

Select Series cushions are one of ROHO's most popular product lines because of the revolutionary ISOFLO Memory Control, which offers shape fitting capabilities with the push of a knob. The ISOFLO Memory Control offers a seated individual an adjustable, quick and simple positioning and stability solution.

"I love my QUADTRO SELECT Cushion and think it is pretty cool! It helps me sit in my wheelchair and protects my skin."

Abel Rose – Waxahachie, TX



FEATURES

The perfect sitting position is only a push away

ROHO's ISOFLO Memory Control, located at the front of all SELECT Series cushions, allows air to flow through the cushion in two ways:



With the ISOFLO unlocked, air can flow freely throughout the entire cushion, similar to a ROHO Single Compartment Cushion.



With the ISOFLO locked, air is contained in four separate quadrants offering stability and positioning capabilities.

PRODUCTS



ROHO QUADTRO SELECT® HIGH PROFILE® Cushion

Approximate Height: 4.25 inches, 10.5 cm

Approximate Immersion: 3.75 inches, 9 cm

Accommodation of pelvic asymmetry, up to approximately 3 inches, 7.5 cm



ROHO QUADTRO SELECT LOW PROFILE® Cushion

Approximate Height: 2.25 inches, 5.5 cm

Approximate Immersion: 1.75 inches, 4 cm

Accommodation of pelvic asymmetry, approximately up to 1 inch, 2.5 cm



ROHO CONTOUR SELECT® Cushion

Approximate Height (Multiple cell height): 2.25 – 4.25 inches, 5.5 – 10.5 cm

Approximate Immersion: 1.75 inches, 4 cm

Accommodation of pelvic asymmetry, up to approximately 1 inch, 2.5cm



WHY ROHO?

Protection

ROHO is the worldwide leader of seating solutions to prevent and treat pressure injuries including deep tissue injury while providing skin and soft tissue protection.

Proof

ROHO has been in over 90 scientific and clinical studies to reinforce and verify that a ROHO DRY FLOATATION® cushion heals, treats and prevents pressure injuries including deep tissue injury.

Adjustability

ROHO DRY FLOATATION cushions, unlike foam products, provide the ability to constantly adjust and conform to an individual's body, now, and as the individual changes throughout the months and years.

Adaptability

ROHO DRY FLOATATION cushions support the shifts and movements that individuals make throughout each day.

Durability

The high quality of ROHO's proprietary neoprene and our ISO certified manufacturing processes instill confidence that the cushion will continue to perform over years.



ROHO QUADTRO SELECT® CUSHIONS

APPLICATIONS

Based on clinical, scientific or engineering evidence, this product may be suitable for individuals who:

- have a history of skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury), skin grafts, and/or myocutaneous flap surgery
- currently have skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- are at risk for skin/soft tissue breakdown
- have normal or impaired sensation
- require increased stability
- require correction of pelvic asymmetry (medial/lateral and/or anterior-posterior)
- require different positioning of the left and right femur

SPECIFICATIONS

INCLUDES: QUADTRO SELECT HIGH PROFILE, MID PROFILE or LOW PROFILE Cushion, two-way stretch cover, hand inflation pump, repair kit, and operation manual.

CONSTRUCTION: Flame-resistant, neoprene rubber (made without latex).

SIZES*: Fits chair width

HIGH PROFILE - 11" to 25"/28 cm to 63.5 cm

MID PROFILE - 13" to 19"/33 cm to 48 cm

LOW PROFILE - 10" to 24"/25.5 cm to 61 cm

*Customization and custom sizes available

WEIGHT: Varies by size, approximately 3.5 lbs./ 1.6 kg

WEIGHT LIMIT: Unlimited (Cushion must be properly sized to the individual)

CUSHION WARRANTY: 36-Month Limited Warranty

COVER WARRANTY: 6-Month Limited Warranty

ROHO CONTOUR SELECT® CUSHION

APPLICATIONS

Based on clinical, scientific or engineering evidence, this product may be suitable for both the individuals listed under the QUADTRO SELECT Cushions and for individuals who:

- benefit from a pre-contoured seating surface
- require lower extremity alignment
- require a greater amount of pelvic control – the cushion's pre-contoured surface allows for aggressive positioning
- require accommodation of pelvic asymmetry, up to approximately 1 in./2.5cm

SPECIFICATIONS

INCLUDES: CONTOUR SELECT Cushion, two-way stretch cover, hand inflation pump, repair kit, and operation manual.

CONSTRUCTION: Flame-resistant, neoprene rubber (made without latex)

SIZES*: Fits chair width – 13" to 24"/33 cm to 61 cm

*Customization and custom sizes available

WEIGHT: Varies by size, approximately 3 lbs./1.4kg

WEIGHT LIMIT: Unlimited (Cushion must be properly sized to the individual)

CUSHION WARRANTY: 36-Month Limited Warranty

COVER WARRANTY: 6-Month Limited Warranty

WHY ROHO?



Expertise

ROHO has been manufacturing cushions for over 45 years. There are currently over one million cushions in use, in over 65 countries worldwide.

Leadership

We continue to actively collaborate with industry experts and associations, and receive council from ROHO's elite Scientific Advisory Board consisting of world leaders in pressure injury research.

Passion

We provide unparalleled customer service and product warranty guarantees to ensure that our customers receive the best products and services possible.

©2016 ROHO, Inc. The following are trademarks and registered trademarks of ROHO, Inc:

CONTOUR SELECT®, DRY FLOATATION®, HIGH PROFILE®, ISOFO®,
LOW PROFILE®, MID PROFILE™, QUADTRO SELECT®, ROHO® and shape fitting technology®.

We recommend consulting with a certified rehab professional when ordering and fitting ROHO products.



The Quality Management System
of ROHO, Inc. is certified to
ISO 9001:2008 and ISO 13485:2003.

permobil

ROHO

Permobil Ltd | Office 312, Britannia Office Quarter Waters Meeting, Britannia Way, Bolton, BL2 2HH, England
Tel. (+44) (0)1484 722 888 | info@permobil.co.uk | permobil.co.uk